

# TRANSCUTANEOUS BILI METER COMPETENCY

Verbalizes  
Understanding

Performs  
Accurately

Patient Use:			
1	States frequency of YCB measurements for newborns.		
2	States frequency where TCB measurements are obtained.		
3	Explains procedure to parents.		
4	Removes the meter from the base and turns on power.		
5	Places the measuring probe perpendicular to the infants sternum or forehead and pushes gently until a click is heard.		
6	Reads the value on display.		
7	Turns off TCB meter.		
8	To test another newborn, presses the "RESET" key.		
9	Cleans measuring probe with Sani-wipes or alcohol swab.		
10	Properly positions meter in charger base.		
11	Documents results in newborn chart.		
Daily Operations:			
1	Removes meter from charging base.		
2	Powers meter from OFF if in ON position.		
3	Turns the meter to the ON position by holding down the "RESET" button		
4	Ensures the letters "CHE" appear in the screen, presses the "RESET" button.		
5	Verifies proper light output by pressing the measuring probe firmly onto the checker until it clicks.		
6	Confirms values fall within the parameter listed on the cover.		
7	Returns meter to normal operating mode by powering the meter OFF and back ON.		
8	Returns meter to charging base.		

Name and Signature:

Date:

Preceptor's Name and Signature:

Date:

Reviewed by:

Date:

